

Product Datasheet

Human C19orf63 ELISA Kit orb405527

Description	Human UPF0510 protein INM02(C19orf63) ELISA kit
Reactivity	Human
Range	78 pg/ml-5000 pg/ml
Sensitivity	19.5pg/ml
Tested Applications	ELISA
Target	C19orf63
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm