

Product Datasheet

Human UCN2 ELISA Kit orb405529

| | |
|----------------------------|--|
| Description | Human Urocortin-2(UCN2) ELISA kit |
| Reactivity | Human |
| Range | 6.25 pg/ml-400 pg/ml |
| Sensitivity | 1.56pg/ml |
| Tested Applications | ELISA |
| Target | UCN2 |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |