

Product Datasheet

Human UTS2D ELISA Kit orb405532

| | |
|----------------------------|---|
| Description | Human Urotensin-2B(UTS2D) ELISA kit |
| Reactivity | Human |
| Sensitivity | 0.039 ng/ml |
| Tested Applications | ELISA |
| Target | UTS2D |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |