

Product Datasheet

Human WISP2 ELISA Kit orb405555

Description	Human WNT1-inducible-signaling pathway protein 2(WISP2) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	WISP2
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	076076
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm