

Product Datasheet

Human XBP1 ELISA Kit orb405558

Description	Human X-box-binding protein 1(XBP1) ELISA kit
Reactivity	Human
Range	12.5 pg/mL-800 pg/mL
Sensitivity	3.12 pg/mL
Tested Applications	ELISA
Target	XBP1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P17861
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm