

## Product Datasheet

### Human ZP3 ELISA Kit [Out of Stock] orb405570

<b>Description</b>	Human Zona pellucida sperm-binding protein 3(ZP3) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 pg/ml-2000 pg/ml
<b>Sensitivity</b>	7.81pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	ZP3
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, cell culture supernates and tissue homogenates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul>