

Product Datasheet

Human KIAA1199 ELISA Kit orb405575

Description	Human KIAA1199 ELISA Kit
Reactivity	Human
Range	1 ng/mL-24 ng/mL
Sensitivity	0.4 ng/mL
Tested Applications	ELISA
Target	CEMIP
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8WUJ3
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm