

Product Datasheet

Human IP3 ELISA Kit orb405578

Description	Human inositol 1,4,5,-trisphosphate(IP3) ELISA kit
Reactivity	Human
Range	5 pg/ml-1000 pg/ml
Sensitivity	2.5 pg/ml
Tested Applications	ELISA
Target	IP3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm