

Product Datasheet

Human 17-OHP ELISA Kit orb405583

Description	Human 17-Hydroxyprogesterone,17-OHP ELISA Kit
Reactivity	Human
Range	0.025 ng/ml-10 ng/ml
Sensitivity	0.025 ng/ml
Tested Applications	ELISA
Target	17-OHP
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm