

Product Datasheet

Human 5-HETE ELISA Kit orb405595

Description	Human 5-hydroxyeicosatetraenoic acid(5-HETE)ELISA Kit
Reactivity	Human
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.81pg/ml
Tested Applications	ELISA
Target	5-HETE
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm