

Product Datasheet

Human 5HT ELISA Kit orb405596

Description	Human 5-Hydroxytryptamine,5HT ELISA Kit
Reactivity	Human
Range	4 ng/mL-800 ng/mL
Sensitivity	2 ng/mL
Tested Applications	ELISA
Target	5HT
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm