

Product Datasheet

Human 6ketoPGF1a ELISA Kit orb405601

Description	Human 6-keto-prostaglandin F1a,6-K-PGF1a ELISA Kit
Reactivity	Human
Range	50 pg/ml-800 pg/ml
Sensitivity	50 pg/ml
Tested Applications	ELISA
Target	6ketoPGF1a
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm