

Product Datasheet

Human 8-iso-PGF2a ELISA Kit orb405605

Description	Human 8-iso prostaglandin F2 α ,8-iso-PGF2a ELISA Kit
Reactivity	Human
Range	125 pg/ml-5000 pg/ml
Sensitivity	125 pg/ml
Tested Applications	ELISA
Target	8-iso-PGF2a
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm