

## Product Datasheet

### Human 8-iso-PGF2a ELISA Kit orb405605

<b>Description</b>	Human 8-iso prostaglandin F2 $\alpha$ ,8-iso-PGF2a ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	125 pg/mL-5000 pg/mL
<b>Sensitivity</b>	125 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	8-iso-PGF2a
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 $\mu$ l • Detection Wavelength: 450 nm