

Product Datasheet

Human PINP ELISA Kit orb405610

| | |
|----------------------------|--|
| Description | Human procollagen I N-terminal peptide,PINP ELISA Kit |
| Reactivity | Human |
| Range | 18.75 pg/ml-1200 pg/ml |
| Sensitivity | 4.68 pg/ml. |
| Tested Applications | ELISA |
| Target | PINP |
| Note | For research use only. |
| Expiration Date | Please enquire. |
| Application Notes | <ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm |