

Product Datasheet

Human alpha Synuclein ELISA Kit orb405628

Description	Human synuclein, alpha (non A4 component of amyloid precursor) oligomer (SNCA oligomer) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Tested Applications	ELISA
Target	SNCA
Note	For research use only.
Sample Types	serum, plasma
Uniprot ID	P37840
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm