

Product Datasheet

Human MS4A1 ELISA Kit orb405650

Description	Human B-lymphocyte antigen (CD20) ELISA Kit
Reactivity	Human
Range	3.12 ng/mL-200 ng/mL
Sensitivity	0.78 ng/mL
Tested Applications	ELISA
Target	MS4A1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P11836
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma and tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm