

Product Datasheet

Human BCL6 ELISA Kit orb405655

Description	Human B-cell leukemia/lymphoma 6,Bcl6 ELISA Kit
Reactivity	Human
Range	18.75 pg/mL-1200 pg/mL
Sensitivity	4.68 pg/mL
Tested Applications	ELISA
Target	BCL6
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P41182
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm