

Product Datasheet

Human CD5L ELISA Kit orb405685

Description	Human CD5 Antigen-like(CD5L) ELISA Kit
Reactivity	Human
Range	6.25 ng/ml-400 ng/ml
Sensitivity	1.56 ng/ml
Tested Applications	ELISA
Target	CD5L
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm