

Product Datasheet

Human CD5L ELISA Kit orb405685

Description	Human CD5 Antigen-like(CD5L) ELISA Kit
Reactivity	Human
Range	6.25 ng/mL-400 ng/mL
Sensitivity	1.56 ng/mL
Tested Applications	ELISA
Target	CD5L
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O43866
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm