

---

## Product Datasheet

### Human Coll2-1 ELISA Kit (orb405698)

<b>Description</b>	Human Coll2-1 ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	47 pg/mL-3000 pg/mL
<b>Sensitivity</b>	11.72 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	Coll2-1
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cell culture supernates
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cell culture supernates • Assay Time: 1-5 hours