

Product Datasheet

Human CORIN ELISA Kit orb405699

Description	Human Corin ELISA Kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Sensitivity	0.156 ng/mL
Tested Applications	ELISA
Target	CORIN
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q9Y5Q5
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm