

---

## Product Datasheet

### Human CX3CR1 ELISA Kit (orb405700)

<b>Description</b>	Human CX3C-chemokine receptor 1,CX3CR1 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25 pg/mL-2000 pg/mL
<b>Sensitivity</b>	7.8 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	CX3CR1
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P49238</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm