

Product Datasheet

Human CXCL3 ELISA Kit orb405701

Description	Human C-X-C motif chemokine 3(CXCL3) ELISA kit
Reactivity	Human
Range	9.37 pg/ml-600 pg/ml
Sensitivity	2.34 pg/ml
Tested Applications	ELISA
Target	CXCL3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm