

## Product Datasheet

### Human IKZF1 ELISA Kit [Out of Stock] orb405711

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human DNA-binding protein Ikaros(IKZF1) ELISA kit   |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 31.25 pg/ml-2000 pg/ml  |
| <b>Sensitivity</b>         | 7.8 pg/ml   |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | IKZF1   |
| <b>Note</b>                | For research use only.  |
| <b>Expiration Date</b>     | Please enquire.   |
| <b>Application Notes</b>   | <ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates, cell lysates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul> |