

Product Datasheet

Human LDHD ELISA Kit orb405715

Description	Human D-Lactate Dehydrogenase,D-LDH ELISA Kit
Reactivity	Human
Range	31.25 mIU/mL-2000 mIU/mL
Sensitivity	7.8 mIU/mL
Tested Applications	ELISA
Target	LDHD
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q86WU2
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm