

Product Datasheet

Human STUB1 ELISA Kit [Out of Stock] orb405716

Description	Human E3 ubiquitin-protein ligase CHIP(STUB1) ELISA kit
Reactivity	Human
Range	23.5 pg/ml-1500 pg/ml
Sensitivity	5.8 pg/ml
Tested Applications	ELISA
Target	STUB1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm