

Product Datasheet

Human EBEA antibody (IgG) ELISA Kit orb405726

Description	Human Epstein-Barr virus early antigen (EBEA) antibody (IgG) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Target	EBEA antibody (IgG)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-3 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450nm