

## Product Datasheet

### Human FLT3LG ELISA Kit (orb405734)

<b>Description</b>	Human FMS-like tyrosine kinase 3 ligand,Flt-3L ELISA Kit...
<b>Reactivity</b>	Human
<b>Range</b>	7.8 pg/mL-500 pg/mL
<b>Sensitivity</b>	1.95 pg/mL
<b>Target</b>	FLT3LG
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>P49771</b>
<b>Expiration Date</b>	Please enquire.