

Product Datasheet

Human GPBAR1 ELISA Kit (orb405745)

Description	Human G protein-coupled bile acid receptor 1(GPBAR1) ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	GPBAR1
Note	For research use only.
Sample Types	serum, plasma
Uniprot ID	Q8TDU6
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm