

Product Datasheet

Human ICOSLG ELISA Kit orb405754

Description	Human ICOS ligand(ICOSLG) ELISA kit
Reactivity	Human
Range	93.75 pg/mL-6000 pg/mL
Sensitivity	23.4 pg/mL
Tested Applications	ELISA
Target	ICOSLG
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O75144
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm