

Product Datasheet

Human IGHG3 ELISA Kit orb405756

Description	Human Ig gamma-3 chain C region(IGHG3) ELISA kit
Reactivity	Human
Range	2.7 ng/mL-2000 ng/mL
Sensitivity	2.7 ng/mL
Tested Applications	ELISA
Target	IGHG3
Note	For research use only.
Sample Types	serum, plasma
Uniprot ID	P01860
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm