

Product Datasheet

Human N-MID-OT ELISA Kit [Out of Stock] orb405781

Description	Human N-MID Osteocalcin,N-MID-OT ELISA Kit
Reactivity	Human
Range	0.45 ng/ml-30 ng/ml
Sensitivity	0.11 ng/ml
Tested Applications	ELISA
Target	N-MID-OT
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm