

Product Datasheet

Human NPHS2 ELISA Kit orb405798

Description	Human podocin ELISA Kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078ng/ml
Tested Applications	ELISA
Target	NPHS2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates and cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm