

Product Datasheet

Human P-Selectin/CD62P ELISA Kit orb405804

Description	Human P-Selectin ELISA kit
Reactivity	Human
Range	0.9 ng/ml-60 ng/ml
Sensitivity	0.225 ng/ml.
Tested Applications	ELISA
Target	P-Selectin/CD62P
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma and tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm