

Product Datasheet

Human AKT1 ELISA Kit orb405806

Description	Human RAC-alpha serine/threonine-protein kinase(AKT1) ELISA kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Sensitivity	0.078 ng/mL
Tested Applications	ELISA
Target	AKT1
Note	For research use only.
Sample Types	serum, plasma, cell lysates
Uniprot ID	P31749
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm