

Product Datasheet

Human protein Red ELISA Kit orb405808

Description	Human protein Red (IK) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Sensitivity	6.25 pg/mL
Tested Applications	ELISA
Target	IK
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q13123
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450nm