

Product Datasheet

Human ROCK2 ELISA Kit orb405812

Description	Human Rho-associated coiled-coil containing protein kinase 2(Rock-2) ELISA Kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Sensitivity	7.8pg/ml
Tested Applications	ELISA
Target	ROCK2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm