

Product Datasheet

Human S100A12 ELISA Kit orb405820

Description	Human S100 calcium binding protein A12/Calgranulin-C(S100A12) ELISA Kit
Reactivity	Human
Range	12.5 ng/ml-800 ng/ml
Sensitivity	3.12 ng/ml
Tested Applications	ELISA
Target	S100A12
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, cell culture supernates and cerebrospinal fluid (CSF)• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm