

## Product Datasheet

### Human SALL4 ELISA Kit orb405829

<b>Description</b>	Human Sal-like protein 4(SALL4) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	37.5 pg/mL-2400 pg/mL
<b>Sensitivity</b>	9.37 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	SALL4
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>Q9UJQ4</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm