

Product Datasheet

Human SLIT3 ELISA Kit orb405833

Description	Human Slit homolog 3 protein(SLIT3) ELISA kit
Reactivity	Human
Range	0.156 ng/ml-10 ng/ml
Sensitivity	0.039 ng/ml
Tested Applications	ELISA
Target	SLIT3
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm