

## Product Datasheet

### Human SMOC2 ELISA Kit (orb405841)

<b>Description</b>	Human SPARC-related modular calcium-binding protein 2(SMOC2) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Sensitivity</b>	0.039 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	SMOC2
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cell culture supernates, tissue homogenates
<b>Uniprot ID</b>	<a href="#">Q9H3U7</a>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm