

Product Datasheet

Human Syndecan 1 ELISA Kit orb405845

Description	Human Syndecan-1/CD138(SDC1) ELISA Kit
Reactivity	Human
Range	4.7 ng/mL-300 ng/mL
Sensitivity	1.17 ng/mL
Tested Applications	ELISA
Target	SDC1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P18827
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm