

Product Datasheet

Human Tau ELISA Kit orb405849

Description	Human Tau proteins ELISA kit
Reactivity	Human
Range	15.6 pg/ml-1000 pg/ml
Sensitivity	3.9 pg/ml.
Tested Applications	ELISA
Target	Tau
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cerebrospinal fluid (CSF)• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm