

## Product Datasheet

### Human TLR2 ELISA Kit orb405854

<b>Description</b>	human toll-like receptor 2,TLR2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312 ng/mL-20 ng/mL
<b>Sensitivity</b>	0.078 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	TLR2
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>O60603</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm