

Product Datasheet

Human alpha 1 Antitrypsin ELISA Kit (orb405880)

Description	Human α 1-Antitrypsin, α 1-AT ELISA Kit
Reactivity	Human
Range	78 nIU/mL-5000 nIU/mL
Sensitivity	19.5 nIU/mL
Tested Applications	ELISA
Target	SERPINA1
Note	For research use only.
Sample Types	serum, plasma
Uniprot ID	P01009
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma • Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm