
Product Datasheet

Human α -MSH ELISA Kit (orb405896)

Description	Human α -melanocyte stimulating hormone(α -MSH) ELISA Kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.039 ng/mL
Tested Applications	ELISA
Target	α -MSH
Note	For research use only.
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm