

Product Datasheet

Human ENO1 antibody (IgA) ELISA Kit orb405900

Description	Human alpha-enolase (ENO1) antibody (IgA) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Target	ENO1 antibody (IgA)
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma• Assay Time: 1-3 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450nm