

Product Datasheet

Human ENO1 ELISA Kit orb405903

Description	Human Alpha-enolase (ENO1/ENO1L1/MBPB1/MPB1) ELISA kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078 ng/ml
Tested Applications	ELISA
Target	ENO1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, tissue homogenates, plasma• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm