

Product Datasheet

Human GLB1 ELISA Kit (orb405913)

Description	Human β -galactosidase, β GAL ELISA Kit
Reactivity	Human
Range	7.8 μ U/mL-500 μ U/mL
Sensitivity	1.95 μ U/mL
Tested Applications	ELISA
Target	GLB1
Note	For research use only.
Sample Types	serum, plasma, urine, tissue homogenates
Uniprot ID	P16278
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, urine, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 μ l • Detection Wavelength: 450 nm